

Cyto Accessories

Cells – fundamental to life

Cytology makes a valuable contribution to medical diagnostics. High cell yield and good cell presentation are essential for the correct evaluation of cytological preparations.

The Hettich cyto system provides cytologists with a wide selection of accessories for preparing slides for the cytological examination of various body fluids.



two-tier adapter **5280** illustrated with slide carriers with fastening ring **1662** and cyto chambers **1663**

The assembly of a complete cyto insert is really easy



suspension 1660 for cyto chambers


suspension 1680 for angle chambers



filter cards



one-funnel chambers



angle chambers



multi-funnel chambers





slide carrier **1662**
slide carrier for 2 chambers **1670**
slide carrier with fastening slider **1470**

DJB Labcare - The UK's Leading Centrifuge Specialist
info@djblabcare.co.uk www.djblabcare.co.uk 01908 612598

Cyto Chambers




cyto chambers for slide carrier 1470

fitting into suspension 1660 and two-tier adapter 5280

cyto chambers		
capacity in ml	8	2 x 8
area in mm ²	240	240
diameter in mm	17.5	17.5
Cat. No.	1471	1475

cyto chambers for slide carrier 1662 and slide carrier for 2 chambers 1670

fitting into suspension 1660 and two-tier adapter 5280

cyto chambers						
capacity in ml	1	2	4	8	3 x 2	4 x 1
area in mm ²	30	60	120	240	3 x 60	4 x 30
diameter in mm	6.2	8.7	12.4	17.5	3 x 8.7	4 x 6.2
Cat. No.	1663	1664	1665	1666	1667	1668
filter cards						
Cat. No.	1675			1676	1677	1678

Angle Chambers

angle chambers for slide carrier 1662

fitting into suspension 1680

angle chambers			
area in mm ²	30	60	120
diameter in mm	6.2	8.7	12.4
Cat. No.	1671	1672	1673
filter cards			
Cat. No.	1696	1697	1698

You can find out how to sediment **liquor cerebrospinalis, urine, bronchial fluids, pleural fluid** and **ascites** on slides, and which accessories to use for obtaining best preparations by consulting Hettich's cyto application notes. These notes can be ordered free of charge or downloaded from our homepage www.hettichlab.com.

slide carriers



slide carrier with fastening slider
Cat. No. 1470



slide carrier with fastening ring
Cat. No. 1662



slide carrier for 2 chambers
Cat. No. 1670

suspensions



lid
Cat. No. 1661



suspension
Cat. No. 1660



suspension
Cat. No. 1680



suspension
Cat. No. 1452

two-tier adapter



Cat. No. 5280

Use of cyto chambers

slide carriers		suspension/adapter			rotor		centrifuge
Cat. No.	qty.	Cat. No.	lid	qty.	Cat. No.	places	
1470 1662 1670	4	1660	1661	4	1624	4	ROTOFIX 32 A UNIVERSAL 320/320 R
	6	1660	1661	6	1626	6	
	4	1452	-	4	1494	4	UNIVERSAL 320/320 R
	4	1452	-	4	1724	4	ROTINA 38/38 R
	6	1660	1661	6	1726	6	
	8	1660	1661	8	1748	8	
	8	5091 / 5280	5095	4	1798	4	
	8	5092 / 5280 / 2 x 5220	5095	4	5624	4	ROTANTA 460/460 R ROTANTA 460 S/460 RS
	8	5091 / 5280	5095	4	5694	4	
	12	5091 / 5280	5095	6	4446	6	
	12	4396/5280	-	6	4256	6	ROTIKA 50 S/50 RS
	12	5091 / 5280	5095	6	4296	6	

Use of angle chambers

slide carriers		suspension/adapter			rotor		centrifuge
Cat. No.	qty.	Cat. No.	lid	qty.	Cat. No.	places	
1662	4	1680	incl.	4	1048	4	MIKRO 22 R
	4			4	1624	4	ROTOFIX 32 A UNIVERSAL 320/320 R
	6			6	1626	6	
	8			8	1648	8	
	6			6	1726	6	ROTINA 38/38 R
	8			8	1748	8	

Disposable Pasteur pipettes for removing the supernatant are available under Cat. No. **0552**.

Please note that RCF values in excess of 1,100 may lead to glass breakage.